

PATENT APPLICATION  
DOCKET NO.: INL-059 (4643/104)

**AQUEOUS SOLUTIONS FOR REDUCING THE RATE OF OXYGEN LOSS, AND  
METHODS THEREOF**

RELATED APPLICATIONS

This application claims the benefit of U.S. Provisional Application Serial Number  
60/398,661, filed on July 26, 2002, the entire disclosure of which is hereby incorporated by  
5 reference.

FIELD OF THE INVENTION

The present invention is related to the field of electrochemical sensors, and calibration  
solutions, particularly to the increased accuracy and effective life span of oxygen sensors.

BACKGROUND OF THE INVENTION

10 In a variety of clinical situations it is important to measure certain chemical  
characteristics of the patient's blood such as pH, hematocrit, the concentration of calcium  
ions, potassium ions, chloride ions, sodium ions, glucose, lactate, creatinine, urea, the partial  
pressure of oxygen (O<sub>2</sub>), carbon dioxide (CO<sub>2</sub>), and the like. These situations range from a  
routine visit of a patient in a physician's office to monitoring of a patient during open-heart  
15 surgery. The required speed, accuracy, and other performance characteristics vary with each  
situation.

Typically, electrochemical sensor systems which provide blood chemistry analysis are  
stand-alone machines or are adapted to be connected to an extracorporeal shunt or an ex vivo  
blood source such as a heart/lung machine used to sustain a patient during surgery. Thus, for  
20 example, small test samples of ex vivo blood can be diverted off-line from either the venous

or arterial flow lines of a heart/lung machine directly to a chamber exposed to a bank of micro-electrodes which generate electrical signals proportional to chemical characteristics of the real time flowing blood sample.

Electrochemical sensor systems are analytical tools combining a chemical or  
5 biochemical recognition component (e.g., an enzyme) with a physical transducer such as a platinum electrode. The chemical or biochemical recognition component is capable of selectively interacting with an analyte of interest and of generating, directly or indirectly, an electrical signal through the transducer. Electrochemical sensor systems play an increasing role in solving analytical and clinical problems, and find applications in the field of medical  
10 diagnostics.

The selectivity of certain biochemical recognition components makes it possible to develop electrochemical sensors that can accurately detect certain biological analytes even in a complex analyte mixture such as whole blood. Despite the high degree of selectivity of electrochemical sensors, the accuracy of such sensors depends on calibrating the sensor  
15 against a calibration solution that includes a known concentration of analyte. If the actual concentration of the analyte in the calibration solution is different from the concentration assumed to be in the solution, then the sensor may be improperly calibrated and the sensor readings may be inaccurate. The oxygen in a calibration solution may be removed by, for example, the oxidation of other components of the solution or contaminants in the solution.  
20 The amount of oxidation that occurs in the calibration solution is unpredictable and, over a period of time, may lead to an oxygen content in the calibration solution different from the oxygen content estimated to exist from the preparation of the solution. A calibration solution

in which the oxygen content is unknown and may be different from the estimated value may lead to an incorrect calibration of the sensor and inaccurate sensor readings.

### SUMMARY OF THE INVENTION

5           One objective of the present invention is to provide a system and a method for increasing the accuracy and effective life span of an electrochemical sensor system, in particular a calibrating solution required for the calibration of an oxygen sensor of the electrochemical sensor system. According to the invention described herein, the addition of choline to an oxygen-containing solution such as a calibrating solution, for example, reduces  
10   the rate at which oxygen is lost from the solution. The reduction in the rate of loss of oxygen content of the solution stabilizes the oxygen content in the solution over a period of time. The stabilized oxygen content in a calibration solution increases the effective life span of the calibration solution's ability to accurately calibrate an oxygen sensor.

          In one aspect of the present invention, a solution for the calibration of an oxygen  
15   sensor includes a selected concentration of choline and a known oxygen content. The concentration of choline is sufficient to reduce a rate of loss of oxygen content in the solution, and the known oxygen content is at a level sufficient to calibrate an oxygen sensor. In one embodiment, the rate of loss of oxygen in the solution is in the range of about 0.05 millimeters of mercury/month (mmHg/month) to about 5.0 mmHg/month. In one  
20   embodiment, the choline in the solution includes choline chloride. In alternative embodiments, other choline salts such as choline hydrogen citrate, choline bitartrate, choline bicarbonate, tricholine citrate, choline ascorbate, choline borate, choline gluconate, choline

phosphate, choline di(choline)sulphate and dicholine mucate may be sources of choline in the solution. In one embodiment, the choline concentration in the solution is in the range of about 5 millimoles/liter (mmol/L) to about 100 mmol/L. In another embodiment, the choline concentration in the solution is about 20 mmol/L.

5           In one embodiment of the invention, the oxygen content in the solution is in the range of about 10 mmHg to about 300 mmHg. In other embodiments of the invention, the oxygen content in the solution is about 100 mmHg or about 180 mmHg.

          In one embodiment of the present invention, the solution for the calibration of an oxygen sensor further includes CO<sub>2</sub>. The solution may also include helium gas. In another  
10   embodiment the solution may include Na<sup>+</sup>, K<sup>+</sup>, Ca<sup>++</sup>, and/or HCO<sub>3</sub><sup>-</sup>. Such compounds may function as calibration standards for electrochemical sensors specific to the compound. In another embodiment, the solution further includes a surfactant. The solution may also include an inert preservative and/or a biological buffer.

          In one aspect of the present invention, a container of calibration solution for  
15   calibrating an electrochemical sensor includes a selected concentration of choline and a known concentration of oxygen. The container for holding the calibration solution is substantially gas-impermeable. According to this aspect of the invention, the concentration of choline is sufficient to reduce a rate of loss of oxygen content in the solution, and the oxygen content in the solution is at a level sufficient to calibrate an oxygen sensor. In one  
20   embodiment, the container includes at least one flexible wall that allows for contraction of a volume of the container. The flexible wall allows the volume of the container to decrease as

the solution is drawn out of the container. In another embodiment, the container enclosing the calibration solution is sealed to prevent a headspace that includes a gas.

In one embodiment of the container, the choline in the solution includes choline chloride. In alternative embodiments of the container, other choline salts such as choline  
5 hydrogen citrate, choline bitartrate, choline bicarbonate, tricholine citrate, choline ascorbate, choline borate, choline gluconate, choline phosphate, choline di(choline)sulphate and dicholine mucate, may be included in the solution. In one embodiment of the container, the choline concentration in the solution is in the range of about 5 mmol/L to 100 mmol/L. In another embodiment, the choline concentration in the solution is about 20 mmol/L.

10 In another embodiment of the container, the oxygen content in the solution is in the range of about 10 mmHg to 300 mmHg. In other embodiments of the invention, the oxygen content in the solution is about 100 mmHg and about 180 mmHg.

In yet another embodiment of the container, the solution for the calibration of an oxygen sensor further includes CO<sub>2</sub>. In other embodiments of the container the solution may  
15 further include helium gas, Na<sup>+</sup>, K<sup>+</sup>, Ca<sup>++</sup>, and/or HCO<sub>3</sub><sup>-</sup>. Such compounds may function as calibration standards for electrochemical sensors specific to the compound. In another embodiment, the solution further includes a surfactant. The solution may also include an inert preservative and/or a biological buffer.

In another aspect, the invention relates to a method of reducing a rate of loss of  
20 oxygen content in a solution. The method includes providing a solution that includes oxygen dissolved in the solution. The method also includes adding choline to the solution, wherein the choline is added in an amount sufficient to reduce the rate of loss of oxygen content in a

solution. In one embodiment, 20 mmol/L choline is added to the solution. In another embodiment, choline is added to the solution in the range of about 5 mmol/L to 100 mmol/L. In another embodiment, choline further includes choline chloride. In yet another embodiment, oxygen in the solution is about 100 mmHg or about 180 mmHg. In another  
5 embodiment, oxygen in the solution is in the range of about 10 mmHg to about 300 mmHg.

In another embodiment of the method, the solution includes a calibration solution. In other embodiments, choline salts such as choline hydrogen citrate, choline bitartrate, choline bicarbonate, tricholine citrate, choline ascorbate, choline borate, choline gluconate, choline phosphate, choline di(choline)sulphate and dicholine mucate may be included in the solution.  
10 In another embodiment, the solution further includes CO<sub>2</sub>. In another embodiment, the solution further includes helium gas. In other embodiments of the method, the solution may further include Na<sup>+</sup>, K<sup>+</sup>, Ca<sup>++</sup>, HCO<sub>3</sub><sup>-</sup>, a surfactant, an inert preservative, and/or a biological buffer.

These and other objects, along with advantages and features of the present invention  
15 herein disclosed, will become apparent through reference to the following description, the accompanying drawings, and the claims. Furthermore, it is to be understood that the features of the various embodiments described herein are not mutually exclusive and can exist in various combinations and permutations.

## 20 BRIEF DESCRIPTION OF THE DRAWING

The foregoing and other objects, features and advantages of the present invention disclosed herein, as well as the invention itself, will be more fully understood from the

following description of preferred embodiments and claims, when read together with the accompanying drawings. The drawings are not necessarily to scale, emphasis instead generally being placed upon illustrating the principles of the invention.

FIG. 1 is a schematic diagram of the components of an electrochemical sensor apparatus including a sensor cartridge with a bank of sensors and a thermal block for accelerated hydration and calibration of the sensors.

FIG. 2 depicts a bar graph of rates of loss of oxygen concentration of solutions in the presence or absence of choline chloride.

10

#### DETAILED DESCRIPTION OF THE INVENTION

The present invention relates to electrodes and electrochemical sensor systems for measuring oxygen levels of aqueous samples including, but not limited to, blood, serum or other body fluids. Specifically, the invention is directed to improving the effective life span of calibration solutions for calibrating an oxygen sensor by stabilizing the oxygen content of the calibration solution over an extended period of time.

15

#### Definitions

In order to more clearly and concisely point out and describe the subject matter which applicant regards as the invention, the following definitions are provided for certain terms used in the following description and claims.

20

As used herein, the term "electrode" refers to a component of an electrochemical device that makes the interface between the external electrical conductor and the internal

ionic medium. The internal ionic medium is typically an aqueous solution with dissolved salts. The medium may also include proteins in a stabilizing matrix.

Electrodes are of three types: working or indicator electrodes, reference electrodes, and counter electrodes. A working or indicator electrode measures a specific chemical species, such as an ion. When electrical potentials are measured by a working electrode, the method is termed potentiometry. All ion-selective electrodes operate by potentiometry. When current is measured by a working electrode, the method is termed amperometry. Oxygen measurement is performed by amperometry. Working electrodes may also function by including an enzyme as part of an enzyme layer that is part of a composite layer that is in close contact with the electrode. The enzyme, which is specific to a particular analyte, may produce hydrogen peroxide, a by-product of the catalytic reaction of the enzyme on the analyte. Hydrogen peroxide is detected by the electrode and converted to an electrical signal. A reference electrode serves as an electrical reference point in an electrochemical device against which electrical potentials are measured and controlled. In one embodiment, silver-silver nitrate forms the reference electrodes. Other types of reference electrodes are mercury-mercurous chloride-potassium chloride or silver-silver chloride-potassium chloride. A counter electrode may act as a sink for the current path.

As used herein, the term “sensor” is a device that responds to variations in the concentration of a given chemical species, such as glucose or oxygen, in a sample, such as a body fluid sample. An electrochemical sensor is a sensor that operates based on an electrochemical principle and requires at least two electrodes. For ion-selective measurements, the two electrodes include an ion-selective electrode and a reference



electrode. Amperometric enzyme electrodes additionally require a third electrode, a counter electrode. Moreover, enzyme sensors based on two electrodes, a working and reference electrode, are also common.

As used herein, the term “ion selective electrode” generally refers to a silver wire  
5 coated with silver chloride in contact with a buffer solution containing a chloride concentration (the inner solution). The buffer solution may be covered with a polymeric ion-selective membrane that is in contact with the test solution. The ion selective membrane typically includes a high molecular weight polyvinyl chloride (PVC), a plasticizer, an ionophore specific to a particular ion, and a borate salt. The surface of the polymeric  
10 membrane is in contact with the test sample on one side and the inner buffer solution on the other side of the membrane.

As used herein, the term “calibration” refers to the process by which the response characteristics of a sensor to a specific analyte are determined quantitatively. To calibrate a sensor, the sensor may be exposed to at least two reagent samples, each reagent sample  
15 having a different, known concentration of an analyte. The responses, i.e., signals, measured by the sensor, relative to the concentrations of the analyte in the two different reagent samples, serve as reference points for measurements of the analyte in samples having unknown concentrations of the analyte.

As used herein, the term “oxygen content” refers to the amount of oxygen dissolved  
20 in a solution. The oxygen content of a solution may be measured in terms of millimeters of mercury (mmHg). The oxygen content of a solution is analogous to the concentration of an analyte in the solution.

As used herein, the term “starting time point” refers to a time point during the life span of an oxygen-containing solution such as, for example, a calibration solution. The starting time point of the solution may, for example, coincide with or immediately follow the solution being placed into a sealed container. The age of the solution is measured from the starting time point of the solution.

As used herein, the term “original oxygen content” refers to the oxygen content of a solution at the starting time point of the solution.

As used herein, the term “rate of loss of oxygen content” refers to a rate that oxygen is lost from a particular solution in a particular container or environment. The rate of loss of oxygen content is experimentally determined by measuring the loss of oxygen content from a solution over a series of time points. Time points may be measured, for example, daily, weekly, or monthly.

As used herein, the term “estimated rate of loss of oxygen content” refers to, for example, a median or average of rates of loss of oxygen concentration from a plurality of batches of a solution. Each batch of the plurality of batches must have identical compositions, identical methods of preparation, and exist in identical containers or environments for the estimated rate of loss of oxygen content calculated from the plurality of batches to be representative of the rate of loss of oxygen content for any one of the batches. Any variation in the compositions, methods of preparation, or containers or environments between the batches may increase the variability of the rates of loss of oxygen content among the batches. An estimated rate of loss of oxygen content calculated from a plurality of non-

identical batches may not be representative of any single batch and thus may not accurately represent the rate of loss of oxygen content for any one solution.

As used herein, the term “estimated oxygen content” refers to a value of oxygen content that is calculated by the electrochemical sensor system in order to, for example, estimate the amount of oxygen in the calibration solution. The estimated oxygen content of the calibration solution is calculated by determining the product of the estimated rate of loss of oxygen content for the calibration solution and the time elapsed since the starting time point of the calibration solution. The resulting product is subtracted from the original oxygen content to determine the estimated oxygen content. The estimated oxygen content is the oxygen content value assigned by the electrochemical sensor system to the calibration solution. The estimated oxygen content of the calibration solution is the value used by the electrochemical sensor system in calibrating the oxygen sensor.

#### Electrochemical Sensor System

Referring to FIG. 1, the exemplary electrochemical sensor system 8 employs a sensor assembly, generally indicated at 10, incorporating a plurality of electrodes 9 adapted to make electrical measurements on a sample, such as a blood sample, introduced to the sensor assembly 10. Blood samples may be obtained by, for example, phlebotomy or are derived on a periodic basis from an extracorporeal blood flow circuit connected to a patient during, for example, open heart surgery. An electrochemical system of a similar type is described in U.S.S.N. 09/872,240, U.S.S.N. 09/871,885, and U.S.S.N. 09/872,247, the entirety of each of the three applications being incorporated by reference herein. An electrochemical system of

a similar type is the GEM Premier 3000 system manufactured by Instrumentation Laboratory, Lexington, Massachusetts 02421.

The electrochemical sensor system 8 includes a disposable cartridge 37. A cartridge of a similar type is set forth in detail in U.S. Patent No. 4,734,184, the entirety of which is  
5 incorporated by reference herein. Electrochemical system 8 also includes a pump 26 for drawing fluids to the electrodes 9, a heater block assembly 39 of a suitably adapted blood chemistry analysis machine, and a microprocessor 40.

In one embodiment of the invention, the electrochemical sensor system 8 incorporates in the cartridge 37 at least two prepackaged containers 14 and 16, each containing a  
10 calibration solution having known values of the parameters to be measured by electrochemical sensor system 8. For purposes of reference, the solution contained within the prepackaged container 14 will be termed calibration solution A, and the solution contained within the prepackaged container 16 will be termed calibration solution B. In another  
embodiment of the invention, the electrochemical sensor system 8 may include a third  
15 prepackaged container 23 containing calibration solution C. Each of the prepackaged containers 14, 16, and 23 contains a sufficient quantity of its respective calibration solution to allow electrochemical sensor system 8 to be calibrated a substantial number of times before the prepackaged container 16 becomes empty, e.g., 1000 times with calibration solution B or 100 times with calibration solution A.

#### 20 Calibration Solutions and Containers

In one embodiment of the invention, a composition of calibration solution A, prepared at 37°C and at atmospheric pressure tonometered with 9% CO<sub>2</sub>, 14% O<sub>2</sub> and 77% helium gas,

is as follows: pH = 6.9; partial pressure of CO<sub>2</sub> (pCO<sub>2</sub>) = 63 mmHg; pO<sub>2</sub> = 100 mmHg; Na<sup>+</sup> = 100 mmol/L; K<sup>+</sup> = 7 mmol/L; Ca<sup>++</sup> = 2.5 mmol/L; HCO<sub>3</sub><sup>-</sup> = 11.5 mmol/L; glucose = 150 milligram/deciliter (mg/dL); lactate = 4 mmol/L; creatine = 0.4 mmol/L; creatinine = 0.3 mmol/L; a biological buffer, e.g., 3-Morpholinopropanesulfonic acid (MOPS) = 100 mmol/L;  
5 ; a surfactant, e.g., polyoxyethylene 23 lauryl ether (Brij<sup>®</sup> 35) = 50 milligram/liter (mg/L); and an inert preservative, e.g., 2-methyl-4-isothiazolin-3-one hydrochloride (MIT) = 0.3 gram/liter (g/L).

In another embodiment of the invention, a composition of calibration solution B, prepared at 37°C and at 700 mmHg absolute pressure tonometered with 27% O<sub>2</sub>, 5% CO<sub>2</sub>,  
10 and 68% helium gas, is as follows: pH = 7.4; pCO<sub>2</sub> = 34 mmHg; pO<sub>2</sub> = 180 mmHg; Na<sup>+</sup> = 140 mmol/L; K<sup>+</sup> = 3.5 mmol/L; Ca<sup>++</sup> = 1.0 mmol/L; HCO<sub>3</sub><sup>-</sup> = 20.0 mmol/L; choline chloride = 20 mmol/L; a biological buffer, e.g., MOPS = 100 mmol/L; a surfactant, e.g., Brij<sup>®</sup> 35 = 50 mg/L; and an inert preservative, e.g., MIT = 0.3 g/L.

In yet another embodiment of the invention, calibration solution C contains an  
15 aqueous solution of Na<sup>+</sup>, K<sup>+</sup>, Ca<sup>++</sup>, and HCO<sub>3</sub><sup>-</sup> salts with concentrations of Na<sup>+</sup> = 140 mmol/L, K<sup>+</sup> = 4 mmol/L, Ca<sup>++</sup> = 0.2 mmol/L, and HCO<sub>3</sub><sup>-</sup> = 13 mmol/L; 15 mmol/L of *m*-phenylenediamine, 20 mmol/L of sulfite; 50 mg/L of a surfactant, e.g., Brij<sup>®</sup> 35; 0.3 g/L of an inert preservative, e.g., MIT; 50 mmol/L of a biological buffer, e.g., MOPS; pH = 7.2; and pCO<sub>2</sub> = 30 mmHg.

20 The reference solution may contain AgNO<sub>3</sub> = 0.5 mmol/L; KNO<sub>3</sub> = 1 mol/L; and a surfactant, e.g., Brij<sup>®</sup> 35 = 150 g/L.

The compositions of the calibration solutions A and B are chosen so that for each of the characteristics measured by electrochemical sensor system 8, e.g., a concentration of  $\text{Na}^+$ , a pair of values are obtained that are spaced over the range of permissible values that are measured by electrochemical sensor system, providing a balanced 2-point calibration for the instrument. The composition of the calibration solution C is selected to enable low level oxygen calibration.

The calibration solutions A and B compositions are prepared by premixing all of the constituents in a certain order, starting with the biological buffer, and then adding the salts of the  $\text{Na}^+$ ,  $\text{K}^+$ ,  $\text{Ca}^{++}$ , and  $\text{HCO}_3^-$  ions, e.g., sodium bicarbonate salt ( $\text{NaHCO}_3$ ). Glucose, lactate, creatine, and creatinine may be added to calibration solution A before or after the addition of the salts. The solution is then tonometered with oxygen and  $\text{CO}_2$  mixed with helium to produce the desired level of  $\text{pCO}_2$  and  $\text{pO}_2$ .

Choline, in the form of, e.g., choline chloride, is added to calibration solution B during the addition of the salts. The concentration of choline may be selected from the range of about 5 mmol/L to about 100 mmol/L. The choline concentration may be, for example, about 20 mmol/L. Choline may also be added in the form of, e.g., one or more of the following compounds: choline hydrogen citrate, choline bitartrate, choline bicarbonate, tricholine citrate, choline ascorbate, choline borate, choline gluconate, choline phosphate, choline di(choline)sulphate and dicholine mucate

The choline in calibration solution B reduces the rate of loss of oxygen content in calibration solution B. The mechanism by which choline reduces the rate of loss of oxygen content in a solution may be as follows. Choline may function by preventing the oxidation of

other compounds present in the solution, and by doing so prevent consumption of the oxygen in the solution in the oxidation reaction. Choline may function to block the oxidation of the compounds by, for example, interacting with the oxygen molecule, the compound to be oxidized, or both. Choline in the solution may also alter the aqueous environment of the solution making it less conducive to an oxidation reaction.

The surfactant and the inert preservative are also added to the solution during the preparation of calibration solution A or B. These two constituents may be added to calibration solution A or B at any point during the preparation of these solutions.

The calibration solution C is prepared by a slightly different procedure. The biological buffer and certain salts (salts of  $\text{Na}^+$ ,  $\text{K}^+$ , and  $\text{Ca}^{++}$  ions) are added to water and the solution is tonometered with helium to bring the  $\text{pO}_2$  to less than 30 mmHg. Then, the remaining chemicals (sulfite, e.g., sodium sulfite; salt of  $\text{HCO}_3^-$  ions, e.g., sodium bicarbonate; and m-phenylenediamine) are added to the solution and the final mixture is tonometered with mixture of  $\text{pCO}_2$  and helium to produce the desired  $\text{pCO}_2$  level. The surfactant and the inert preservative are also added to the solution during the preparation of calibration solution C. These two constituents may be added to calibration solution C at any point during the preparation of this solution.

The temperature and pressure at which the calibration solutions are prepared and their method of packaging may be selected such as to preclude the possibility of dissolved gases going out of solution into the container and to minimize the tendency for gases to permeate through the extremely impermeable materials of the container. The calibration solutions may be packaged with the solutions completely filling the containers, so that there is no

headspace, by evacuating the containers prior to filling in a manner which will be subsequently described. Elimination of a headspace in the container is important to prevent oxygen dissolved in the calibration solution from diffusing out of solution into the headspace. Diffusion of oxygen from the solution into the headspace may cause the rate of loss of  
5 oxygen concentration in the solution to increase in an unpredictable manner.

By filling the calibration solution into the evacuated and flexible walled prepackaged containers 14, 16, 23 at elevated temperatures and subatmospheric pressure, the solution will not have a tendency at a lower use temperature to outgas and thus produce gas bubbles in the container. Were outgassing to occur, the concentrations of the gases in the solution would be  
10 affected, possibly creating an unpredictable oxygen concentration in the solution. Similarly, the calibration solutions are preferably not packaged at a pressure that is too low, i.e., not below about 625 mm of mercury, because the absorptive capacity of the solution for gases may increase as the packaging pressure decreases, and below that pressure value the absorptive capacity of the solution may be sufficiently high to draw gases in through the  
15 slight inherent permeability of even the most gas impervious flexible packaging material, over long periods of time. Accordingly, a packaging pressure in the range of 625 - 700 mm of mercury is preferred.

In one embodiment, the calibration solution may be prepared at a temperature in excess of its intended use temperature, thereby reducing the tendency for outgassing of the  
20 dissolved gases at the lower temperature. This cooperates with the reduced pressure packaging to minimize the possibility of outgassing. Calibration solutions A and B may each be, for example, prepared at a temperature above their respective intended use temperature at



a controlled pressure close to atmospheric pressure. Through use of elevated temperature (e.g., 37°C), the calibration solutions may be prepared at about atmospheric pressure without any possibility of subsequent microbubbles within the container or gas transfer through the container when packaged in a zero headspace flexible gas impervious container.

5           The calibration solution prepackaged containers 14, 16, 23 may be formed from envelopes. These envelopes may be formed, for example, of rectangular sheets, heatsealed at the edges and heatsealed at one corner to an inlet stem of the valve 18 that is used for filling purposes. In the illustrated preferred embodiment, the prepackaged containers 14, 16, and 23 and the prepackaged container lines 20, 22, and 25 are formed in a unitary cluster with the  
10   valve 18 so that gas phase dead space in the lines 20, 22, 25 is thereby avoided. In a preferred procedure for purging and filling the envelope bags, a bag is first evacuated and then filled with the prepared solution. The bag is then shaken while the excess solution continually flows out of the bag. This process removes any residual gas bubbles from the bag. The solution is then sealed in the container. In a preferred embodiment, prepackaged  
15   containers 14, 16, and 23 include gas-impermeable containers that substantially prevent oxygen from diffusing into or out of the container, either of which would alter the estimated oxygen content used for calibrating the oxygen sensor.

          The calibration solutions in the prepackaged containers 14, 16, and 23 have excellent stability and a long shelf life. When the calibration solutions A, B, and C are at use  
20   temperature and atmospheric pressure, there is minimal possibility of any outgassing from the liquid to form gas bubbles within the prepackaged containers 14, 16, and 23. To further improve the shelf life of the calibration solution B and all embodiments of the invention that

include calibration solution B, the fluctuation in oxygen content in the calibration solution B may be stabilized. Choline chloride may be added to calibration solution B to reduce the rate of loss of oxygen content of the solution, thus stabilizing the oxygen content in the solution.

#### Oxygen Sensor Calibration

5           The oxygen sensor of the electrochemical sensor system may be calibrated at an upper calibration point with reference to the calibration solution B and at a lower calibration point with reference to the calibration solution C. The calibration solution C has an oxygen concentration of approximately zero as a result of the oxygen scavenger sulfite present in the solution. In one embodiment according to the invention, the calibration solution B is  
10 originally formulated to contain an original oxygen content of 180 mmHg.

          The oxygen content of calibration solution B tends to decrease over time and thus differs from the original oxygen content at a time point following the starting time point of the solution. This decrease in the oxygen content may be measured by the rate of loss of oxygen content of calibration solution B. The underlying cause of the rate of loss of oxygen  
15 content of the calibration solution B may be the diffusion of oxygen out of the solution or an oxidation reaction occurring in the solution. The rate of loss of oxygen content of the solution may also vary unpredictably between different batches of an identical calibration solution due to, for example, variations in materials from which the calibration bag is made, variations in the process of sealing the calibration bag, variations in storage temperature  
20 during the shelf life, and impurities that may be present in the calibration solution B.

          An unpredictable rate of loss of oxygen content in the calibration solution B may cause the estimated oxygen content to differ from the oxygen content of calibration solution

B. This difference in values between the estimated oxygen content and the actual oxygen content of calibration solution B may result in the miscalibration of the oxygen sensor.

Calibration of an oxygen sensor occurs by assigning the value of the estimated oxygen content of calibration solution B to the measurement by the oxygen sensor of the oxygen

5 content of calibration solution B. The degree to which the estimated oxygen content of calibration solution B differs from the oxygen content of calibration solution B is the degree of error of the calibration of the oxygen sensor.

Choline, when added to a solution containing oxygen, decreases the rate of loss of oxygen content from a solution, and thus stabilizes the oxygen content in the solution. A  
10 stabilized rate of loss of oxygen content reduces the potential degree of difference between the estimated oxygen content and oxygen content in the calibration solution B, thus allowing for an increased accuracy in the calibration of the oxygen sensor. The increased accuracy of the calibration of the oxygen sensor also extends the useful life span of the calibration solution B. The extended useful life span of calibration solution B is due to the increased  
15 amount of time that must elapse for a difference between the estimated oxygen content and oxygen content of calibration solution B to become sufficiently large to prevent an accurate calibration of the oxygen sensor.

#### EXEMPLIFICATION

Test samples of different calibration solutions labeled sample 1, sample 2, sample 3,  
20 sample 4, and sample 5, were prepared as follows. Samples 1 and 2 were identical to calibration solution B without added choline chloride. Samples 1 and 2 contained 15 and 30 units/ml of heparin, respectively. Heparin is a compound known to increase oxygen decay

rates, and added to samples 1-5 functions to simulate an oxygen degrading compound that may be found in calibration solution B. Samples 3, 4, and 5 were identical to calibration solution B with 20 mmol/L of added choline chloride. Sample 3 and 4 included 30 units/ml of heparin and were otherwise identical in their composition. Sample 5 differed from samples 3 and 4 by including 60 units/ml of heparin. Thus, samples 2, 3, and 4 had identical concentrations of heparin, samples 3 and 4 included choline chloride, and sample 2 did not include choline chloride.

For the test, the oxygen concentrations of the five sample types were analyzed once a month for a period of 6 months. Referring to FIG. 2, a bar graph including the rate of loss of oxygen content of the five test samples is depicted. Each of the five bars depicted in FIG. 2 represents the average of three identical independently prepared samples. A significant reduction in rates of loss of oxygen content is present in samples that included choline chloride. The choline chloride in calibration solution B as part of an electrochemical sensor system solution functions similarly in lowering oxygen decay rates and thus increases the effective life span of the calibration solution B.

Choline added to calibration solution B at a concentration of 20 mmol/L or greater is sufficient to reduce the rate of loss of oxygen content in calibration solution B sufficiently to extend the effective life span of calibration solution B. Concentrations of choline less than 20 mmol/L may not reduce the rate of loss of oxygen content in calibration solution B to the same extent as a concentration of 20 mmol/L. A concentration of choline of 100 mmol/L or lower is sufficient to minimize the potential interference of choline with analytes in calibration solution B, or with electrodes in which calibration solution B makes contact. The

addition of choline to calibration solution B in concentrations greater than 100 mmol/L increases the potential that interference may occur between choline in calibration solution B and analytes in calibration solution B, or between choline in calibration solution B and electrodes which contact calibration solution B.

5           The function of choline in calibration solution B is readily applied to other oxygen-containing solutions. An aspect of the invention relates to a method of stabilizing the oxygen content of a solution by the addition of choline to the solution. This method is applicable to any type of solution that contains dissolved oxygen and it is beneficial for the oxygen content of the solution to remain stable. This includes solutions related to clinical, pharmaceutical,  
10   and industrial applications.

          The invention may be embodied in other specific forms without departing from the spirit or essential characteristics thereof. The foregoing embodiments are therefore to be considered in all respects illustrative rather than limiting of the invention described herein. The scope of the invention is thus indicated by the appended claims rather than by the  
15   foregoing description, and all changes which come within the meaning and range of equivalency of the claims are therefore intended to be embraced herein.

What is claimed is: